## Bissregistraratmasse *Genie*<sup>®</sup> Bite

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#### Einleitung

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**Directions for Use – ENGLISH**

#### Introduction

*Genie*<sup>®</sup> Bite is an addition-cured polyvinyl siloxane (addition silicone) material specifically designed for bite registration procedures. *Genie*<sup>®</sup> Bite is packed exclusively in 50 ml safety cartridges and is available in two (2) set times: Fast (90 seconds) and Super Fast (60 seconds).

*Genie*<sup>®</sup> Bite is a highly thixotropic material that resists slumping and stays where placed. To ensure accurate model re-positioning, both types (set times) of *Genie*<sup>®</sup> Bite set to a high Shore “A” hardness. *Genie*<sup>®</sup> Bite is designed to be easy to trim when set and has a pleasant mint flavor.

#### Indications for Use

*Genie*<sup>®</sup> Bite is indicated use as a dental bite registration material and/or a dental impression material. *Genie*<sup>®</sup> Bite may be used with any bite registration technique, including digital scanning techniques. Additionally *Genie*<sup>®</sup> Bite can be used as a matrix material during the temporary crown and bridge process.

#### Contraindications

None known.

#### Warnings

Do not use as a temporary liner in dentures. Use of material in this way may cause irritation of the oral mucosa.

#### Precautions

- Use only as directed.
- Avoid contamination of *Genie*<sup>®</sup> Bite with latex, vinyl or nitrile gloves which have powder or starch lubricant on the glove surface. Lubricants commonly used on gloves will interfere with the set of *Genie*<sup>®</sup> Bite. Handle retraction cords with cotton pliers and avoid touching prepared teeth with gloved fingers.
- Rinse the area thoroughly when using certain types of hemostatic agents, including Superoxol (hydrogen peroxide), ferric sulfate and aluminum chloride. These can interfere with the setting reaction.
- Do not use in combination with polyether, polysulfide or condensation silicone impression materials (using *Genie*<sup>®</sup> Bite as a bite registration after making an impression with one of the above listed materials is acceptable).
- Before use, make sure *Genie*<sup>®</sup> Bite is at room temperature.
- To ensure optimal working properties the material should be homogeneous (streak free).
- If any question regarding interference with the setting reaction exists, conduct a timed test mix to ensure compatibility.

#### Storage and Shelf Life

*Genie*<sup>®</sup> Bite has a shelf life of 30 months when stored at room temperature (50º- 80°F /10º- 27°C). Do not use after the expiry date marked on the package.

## Adverse Reactions

Do not use as a temporary liner in dentures. Use of material in this way may cause irritation of the oral mucosa. (See Warnings)

## Step-By-Step Instructions For Use

**Bleeding the cartridge is a critical step to ensure the best results.**

- Load the cartridge into the dispenser.
- Remove the color coded cap from the cartridge and discard.
- Squeeze the trigger slowly on the dispenser to extrude the material onto a mixing pad to ensure both base and catalyst are flowing evenly.
- Immediately affix the appropriate mixing tip to the cartridge (6.5 mm green tips are supplied, but 4.2 mm yellow tips will also work).
- Squeeze the trigger on the dispenser, then extrude a small pea-sized amount to ensure proper mixing of base and catalyst.
- Dispense as desired. Do not remove the mixing tip. Allow bite material to harden in the mixing tip to serve as a cap on the cartridge until the next use.
- To perform bite registration, extrude the *Genie*® Bite directly on the occlusal surfaces of the posterior teeth.
- If performing an edentulous recording using a central bearing device, close the patient into centric and dispense *Genie*® Bite along the buccal surfaces of the teeth, ensuring material flows onto the occlusal surfaces of the denture.
- Once the bite registration has been secured and disinfected, it may be trimmed using rotary instruments or a sharp knife.

## Disinfection

- We recommend using our *Sporox*® II sterilant/high-level disinfectant solution to disinfect *Genie*® Bite, however any standard disinfecting solution or spray can be used.
- After removal from the mouth, rinse the bite registration under cold water to remove blood and saliva.
- To disinfect using *Sporox*® II, soak the registration for 30 minutes in the disinfectant solution.
- Rinse the bite registration thoroughly under running water and dry with compressed air.

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