Directions for Use – English
Genie Magic Mix® Ultra Hydrophilic Impression Material
Rapid Set / Standard Set

REF 78800 Genie Magic Mix®, Regular Body Rapid Set, Standard Package
REF 78805 Genie Magic Mix®, Regular Body Standard Set, Standard Package
REF 78810 Genie Magic Mix®, Heavy Body Rapid Set, Standard Package
REF 78815 Genie Magic Mix®, Heavy Body Standard Set, Standard Package
REF 78820 Genie Magic Mix®, Regular Body Rapid Set, Bulk Package
REF 78825 Genie Magic Mix®, Regular Body Standard Set, Bulk Package
REF 78830 Genie Magic Mix®, Heavy Body Rapid Set, Bulk Package
REF 78835 Genie Magic Mix®, Heavy Body, Standard Set, Bulk Package
REF 78840 Genie Magic Mix®, Dynamic Mixing Tips, Refill Package
REF 78845 Genie Magic Mix®, Bayonet Ring for Mixing Tips, Refill Package

Introduction
Genie Magic Mix® is an elastomeric dental impression material based upon addition cured vinyl polysiloxane chemistry. The product is supplied in 380 ml cartridges in 2 viscosities (regular and heavy). Both viscosities are available in both Rapid Set (2:10 total set time) and Standard Set (4:30 total set time). Genie Magic Mix® exhibits excellent hydrophilicity and ideal trixotropic characteristics for any impression technique. The Genie Magic Mix® cartridge will fit into all commercially available 5:1 ratio mixing machines including the 3M/ESPE Pentamix® and the DMG MixStar®.

Indication for Use
Genie Magic Mix® is indicated for use as a dental impression material. The material is intended for use with all dental impression techniques including dual-phase, mono-phase and double impression techniques. Genie Magic Mix® is fully compatible with Genie® Ultra Hydrophilic Impression Material delivered in 50 ml cartridges.

Composition
Genie Ultra Hydrophilic Impression Material is an addition-cured vinyl polysiloxane. It is supplied as 2 pastes in a 5:1 base:catalyst mixing ratio. The product is packed in 380 ml cartridges.

Contraindications
Do not use with polyether, polysulfide impression materials or condensation-type silicone materials. Do not use as a temporary soft liner in dentures.

† Trademarks are the property of their respective owners.

Precautions
Use only as directed. Lotions and lubricants commonly used on hands and gloves will interfere with the set of addition silicones, including Genie Magic Mix®. Use caution to avoid contamination of the Genie Magic Mix® material with latex, vinyl or nitrile gloves which have powder or starch lubricant on the glove surface. Handle retraction cords with cotton pliers and avoid touching prepared teeth with gloved fingers. Certain types of hemostatic agents, including Superoxol, ferric sulfate and aluminum chloride can interfere with the setting reaction. If any question regarding interference with the setting reaction exists, conduct a timed test mix to ensure compatibility.

Pouring the Cast
Genie Magic Mix® is fully compatible with all dental gypsum and epoxy materials. To avoid surface porosity in the cast, delay pouring gypsum for 30 minutes and epoxy materials for 60 minutes. Genie impressions may be copper or silver-plated. The completed impression may be stored for 14 days prior to pouring without loss of accuracy.

Shelf-Life
Genie Ultra Hydrophilic Impression Material has a shelf life of 2.5 years when stored at room temperature (50-80°F / 10-27°C). Do not use after the expiry date marked on the package has passed.

Disinfection
We recommend using our Sporox® sterilant/high-level disinfectant solution to disinfect Genie Magic Mix® impressions, however any standard disinfecting solution or spray can be used. After removal from the mouth, rinse the impression under cold water. To disinfect, soak the impression for 30 minutes in disinfectant solution. Thoroughly rinse and blow dry.

Usage Instructions
1) Remove the outlet cap from the cartridge by gently pulling the retainer tab up and lifting the cap from the cartridge (Figure 1).
2) Load the cartridge into the dispensing machine (Figure 2).
3) Engage machine plunger until it contacts the cartridge pistons (Figure 3).
4) Activate the mixing machine and bleed the cartridge until the base and catalyst are evenly dispensing (Figure 4).
5) Affix the dynamic mixer to the cartridge. Ensure that the machine's drive shaft is engaged in the mixer drive socket. Secure the mixer to the cartridge using the blue bayonet ring (turn clockwise until it locks). Ensure the bayonet ring is fully engaged (Figure 5).
6) Activate the machine and dispense Genie Magic Mix® as desired. (Figure 6).

Genie Magic Mix® is designed to be used with virtually all impression techniques. Genie Magic Mix® is also fully compatible with cartridge-dispensed Genie® Ultra Hydrophilic Impression Material. The selection of the correct viscosity and tray combination is a matter of personal preference. Genie® is fully compatible with all common tray materials. When using a tray without perforations (i.e. a full-custom tray), a tray adhesive is required. We have assembled some of the more common impression scenarios below:

Scenario: Full Arch Impression
Tray: Custom tray, preformed plastic or metal tray
Recommended Materials: Light or regular body wash / regular or heavy body tray
Technique: Either dual-phase or mono-phase.
Details: Genie® is ideal for either a dual-phase or mono-phase technique when making full arch impressions. If using a standard dual-phase technique, either the light body or regular body material may be used as the wash material - the choice is dependent solely on personal preference - with heavy body material used in the tray. Alternatively, regular body material - may be used as both the wash material and tray material in a mono-phase impression. Genie® impression materials work flawlessly together and can be combined in any manner the user desires.

Scenario: Partial Arch Impression
Tray: Custom tray, preformed plastic or metal tray
Recommended Materials: Light or regular body wash / regular or heavy body tray
Technique: Either dual-phase or mono-phase.
Details: If making a partial arch impression, the technique is similar to a full arch impression except for the selection of the tray. As in the case of a full arch impression, partial arch impressions can be made using either dual phase or mono phase techniques.

Scenario: Dual Arch Impression
Tray: Dual arch tray (i.e. a "triple-tray")
Recommended Materials: Light or regular body wash / regular or heavy body tray
Technique: Either dual phase or mono phase.
Details: The dual-arch impression technique is one of most used in dentistry today and the Genie® system is ideal for this popular technique. The dual-arch technique is ideal for single unit crowns and intra-coronal restorations. To ensure rigidity of the tray, we recommend using light or regular body material as the wash and heavy body material in the tray. If the clinician desires, a mono phase technique can be employed using exclusively regular body Genie® material.

Scenario: Double Impression
Tray: Preformed plastic or metal tray
Recommended Materials: Extra light or light body wash, putty tray material.
Technique: The double impression technique is one of the best-established crown and bridge impression techniques. The technique involves using a preformed plastic or metal tray. Genie Putty is kneaded and placed into the tray. The tray is inserted into the patient's mouth and allowed to set. Upon removal, the impression is thoroughly washed and dried. The inter-dental "tags" of Genie material are then removed using a carbide bur. With the same bur, remove approximately 1-2 mm of material around the impression of the prepared teeth and a lesser amount on adjacent teeth. Wash and dry the impression again thoroughly. Isolate the prepared teeth and apply either extra light or light body wash material to the preparations and adjacent teeth. Reseat the altered impression and allow the wash material to set. An alternate usage is a dual-phase putty technique where wash material is injected around the prepared teeth and putty is used in the tray in a single-stage impression (the wash and putty are allowed to cure together).

ULTRA HYDROPHILIC IMPRESSION MATERIAL

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