Directions for Use

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Introduction

Silgimix® is an elastomeric dental impression material based on addition-cure vinyl polysiloxane chemistry. Silgimix® is specifically formulated to replace alginate impression material in the dental practice, giving the user clear benefits of silicone technology which include:

- Long-term dimensional stability – no need to pour immediately
- Multiple pours may be made from the same impression
- Impression can be reliably disinfected
- Impression can be transported long distances
- Impression can be scanned for a variety of CAD/CAM applications.

Silgimix® 380 ml cartridge

Silgimix® is designed especially for taking impressions for study models, orthodontic models, counter impressions, partial denture impressions and temporary crown fabrication. It's unique combination of optimum tear strength, ideal setting hardness, fluffy, non-runny consistency, pleasant mint scent and value pricing make Silgimix® the best choice to replace messy alginate impression materials.

Silgimix® 50 ml cartridge

Silgimix® is supplied in both 50 ml cartridges in a 1:1 mixing ratio and in 380 ml cartridges in a 5:1 mixing ratio. The 380 ml cartridge system features more material and less waste than competitive foil bag systems. The Silgimix® 380 ml cartridge will fit into all commercially available 5:1 ratio mixing machines.

Silgimix® 380 ml cartridge

Indications for Use

Silgimix® is indicated for use as a dental impression material. The material is intended primarily for taking impressions for study models, orthodontic models, counter impressions, partial denture impressions and temporary crown fabrication.

Silgimix® 380 ml cartridge

Technical Specifications

Classification: ISO 4823, Type 2, Medium Consistency
Working Time: 1:00
Intraoral Setting Time: 3:00
Strain-in-Compression: 3-6%
Final Hardness (Shore “A” scale): 45 @ 24h
Linear Dimensional Change: 0.08%
Composition: Ploydimethylsiloxane polymer; polymethylhydrogen Siloxane; Silicon Dioxide; Starch; Calcium Carbonate; Organic Platinum Complex; Surfactant; FD & C Blue #2 Pigment; Peppermint Oil
Color: Light blue (base) and white (catalyst)
Scent: Wild Mint

Contraindications

Silgimix Material is contraindicated for use with patients who have a known hypersensitivity to peppermint oil or any of the components.

Warnings

Do not use Silgimix Material as a temporary reliner. Use of the material in this manner may cause irritation to the oral mucosa. If symptoms occur discontinue use immediately. If irritation does not subside consult a physician.

Precautions

Use only as directed. Hand lotions, skin lubricants, powdered glove and latex may interfere with the setting reaction of addition silicone materials, including Silgimix. Use caution to avoid contamination of the impression material with these items.

Silgimix® 380 ml cartridge

Silgimix material should not be intermixed with, or used in conjunction with, other manufacturers’ vinyly polysiloxanes or any polyether, polysulfide or conventional (condensation-cured) silicones.

Tray Selection

Silgimix® is fully compatible with all commercially available dental impression tray systems including metal trays, plastic trays and partial arch trays. We recommend use of a VPS tray adhesive with all plastic trays.

Disinfection

We recommend using our Sporox® II sterilant/high-level disinfectant solution to disinfect Silgimix® impressions, however any standard disinfecting solution or spray can be used. After removal from the mouth, rinse the impression under cold water. To disinfect, soak the impression for 30 minutes in disinfectant solution. Thoroughly rinse and blow dry.

Pouring the Cast

Silgimix® impressions are fully compatible with all dental gypsum and epoxy materials. To avoid surface porosity in the cast, delay pouring gypsum for 30 minutes and epoxy materials for 60 minutes. The completed impression may be stored for 14 days prior to pouring without loss of accuracy.

Usage Instructions (380 ml cartridge)

- Remove the outlet cap from the cartridge by gently pulling the retainer tab up and lifting the cap from the cartridge (Figure 1).
- Load the cartridge into the dispensing machine (Figure 2).
- Engage machine plunger until it contacts the cartridge pistons (Figure 3)
- Activate the mixing machine and bleed the cartridge until the base and catalyst are evenly dispensing (Figure 4).
- Affix the dynamic mixer to the cartridge. Ensure that the machine’s drive shaft is engaged in the mixer drive socket. Secure the mixer to the cartridge using the blue bayonet ring (turn clockwise until it locks). Ensure the bayonet ring is fully engaged (Figure 5).
- Activate the machine and dispense Silgimix® as desired into impression trays and/or syringes. (Figure 6).

Usage Instructions (50 ml cartridge)

- Load the cartridge into an S-50 dispenser.
- Remove the green colored protective cap from the cartridge and discard.
- Seat the dispenser plunger against the cartridge pistons.
- Squeeze the trigger on the dispenser to bleed a small amount of material from the cartridge, ensuring that both sides are flowing evenly. Discard the bled material.
- Affix a green (6.5 mm) mixing tip onto the cartridge.
- Squeeze the trigger on the dispenser to extrude the impression material. Discard the first, pea-sized amount of material from the mixing tip.
- Dispense as desired.

Silgimix® 50 ml cartridge

Silgimix® is designed to be used anywhere alginate would be used. The additional benefits of addition silicone technology and optimum physical properties make Silgimix the ideal choice for all situations.

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