


**Genie**<sup>®</sup>  
ULTRA HYDROPHILIC IMPRESSION MATERIAL

**MADE TO PERFORM.  
PRICED TO PLEASE.**



**Sultan** Healthcare

# RELIABLE PERFORMANCE— NO MATTER THE TECHNIQUE.

Discover Genie® brand from Sultan Healthcare—a dental materials leader with 150+ years of expertise.

After decades of proven performance, the Genie® brand delivers the impression and bite materials dentists reach for—for detailed, precise impressions that save time and money\*. That's because it's chemically engineered to produce reliable results, every time, no matter the use or technique.


Made to perform. Priced to please.



\*Price Positioning Analysis, Impression and temp materials, Keystone 2024

# A GREAT FIRST IMPRESSION IS PRICELESS.

Engineered for precise dental impressions, Genie® impression material delivers exceptional accuracy, dimensional stability and ease of use. Thanks to its versatility, it ensures a great experience for clinicians and patients alike.



## Why choose Genie® impression material?



### Versatility

Choose from a variety of options:  
 5 viscosities (putty, heavy, regular, light, extra light),  
 2 setting times (rapid or standard),  
 3 delivery systems (machine-mixed, automix-cartridge, hand-mixed putty).



### Flowability

Genie® materials flow readily into the tiny crevices of cavity preparations to accurately capture margin details.\*



### Patient comfort

Genie® comes in a pleasant berry flavor and has a rapid setting time for patient comfort.



### Visual confidence

Color combination and contrast help clinicians feel more confident in their impression interpretations.

\*Material flow "Shark Fin Study" test (26mm) performed by Tufts University. Unpublished study. Data can be requested at [customer.service@sultanhc.com](mailto:customer.service@sultanhc.com)

# HEAR WHAT GENIE® USERS HAVE TO SAY

For decades, Genie® material has proven its worth to clinicians through its high performance and customer satisfaction—all at a great value. But don't take our word for it—see what Genie® users have to say.\*

## **They're brand loyal**


56.7% of Genie® users are loyal to the brand—for 10 years or more!

## **Genie® is a product dentists recommend**

The majority of respondents started using Genie® thanks to a colleague recommendation

## **They appreciate its benefits**

Survey respondents recognized reliability, price-to-performance ratio and versatility as the top three benefits of Genie®.



# WITH GENIE<sup>®</sup>, YOU'VE GOT OPTIONS

Every case—and clinician preference—is unique. That's why Genie<sup>®</sup> comes in a variety of viscosities, setting times and delivery systems.



## Viscosity options

Depending on your impression technique and personal preference, choose from the following:

- Genie<sup>®</sup> flowable viscosities in “Extra light”, “Light” + “Regular” for syringe application around the tooth
- Firmer viscosities in “Putty”, “Heavy” + “Regular” are for use as tray material

## Delivery system options

Genie<sup>®</sup> is available in a range of delivery systems:

- 445 gram jars
- 380 ml cartridges (for use with mixing machines)
- 50 ml cartridges (for use in combination with a dispenser gun)

Material	ISO Type	Viscosity	Use	Color
Putty	0	Putty	Tray	●
Heavy	1	Heavy Body	Tray	●
Regular	2	Medium Body	Tray-Wash	●
Light	3	Light Body	Wash	●
Extra Light	3	Extra Light Body	Wash	●

### Genie<sup>®</sup> Insider Tip

“Regular” viscosity is the only viscosity that can be used for both syringe and tray application—making it the recommended viscosity for simpler cases and the monophasic technique. For more challenging cases, extra light or light are recommended.

# GET A BETTER BITE.

Registering the bite accurately is crucial for the success of a restoration<sup>1</sup>. In addition to the practitioner's skill and technique, the material used can also have a decisive impact on the accuracy of bite registration<sup>2</sup>.

Genie<sup>®</sup> Bite is specifically designed to help dentists achieve accurate bite registration. It starts off with a light, fluffy consistency that stays where you put it without slumping. Once set, it becomes rock-hard, so it captures and keeps great detail without compressing or flexing during model mounting.



## Why choose Genie<sup>®</sup> Bite?



### Hardness

Genie<sup>®</sup> Bite has a quick, rock-hard set to a Shore "A" hardness of 90+. This helps to resist any distortion and maintain an accurate bite.



### Patient comfort

With a pleasant mint flavor and fast (90-second) and super fast (60-second) setting times, patients remain comfortable.



### Ideal consistency

Genie<sup>®</sup> Bite has a light, fluffy consistency that doesn't cause resistance, interference with occlusion or significant thickness between the patient's arches<sup>3</sup>.



### Ease of handling

50 ml automix cartridges deliver ease of handling during preparation and use of material<sup>4</sup>.

<sup>1</sup> Tejo K S., Kumar A G, Kattimani V S, Desai P D, Sandeep N, Chaitanya K. A comparative evaluation of dimensional stability of three types of interocclusal recording materials- an in-vitro multi-centre study'. Head & Face Medicine. 2012; 8(27): 1-9. doi: 10.1186/1746-160X-8-27  
<sup>2</sup> Nagrath R, Lahori M, Kumar V, Gupta V. A comparative study to evaluate the compression resistance of different interocclusal recording materials: an in vitro study. J Indian Prosthodontic Society. 2014; 76-85. doi: 10.1007/s13191-014-0369-8 <sup>3</sup> Mobilio M, Catapano S. Effect of experimental jaw muscle pain on occlusal contacts. Journal of Oral Rehabilitation. 2011; 38(6): 404-409. doi: 10.1111/j.1365-2842.2010.02173. <sup>4</sup> Powers J, Wataha J. Dental Materials foundations and applications. 2017; 15

# YOU SPOKE. WE LISTENED.

According to recent market research, 50% of clinicians surveyed said they strongly prefer an open mix tip system.\* Customer satisfaction is very important to us—so we’re going back to our open mixing system. That means our cartridges now work with the most common mix tips on the market—with the same impression material you know and love.

## Find the correct mix tip

Product	Use	Cartridge	Mix tip
Genie® Putty	Tray		
Genie® Heavy	Tray	50 ml/380 ml**	
Genie® Regular	Tray/Wash	50 ml/380 ml**	
Genie® Light	Wash	50 ml	
Genie® Extra Light	Wash	50 ml	
Genie® Bite	Bite	50 ml	



NEW MIX TIPS  
ALSO AVAILABLE  
FOR 50 ML  
DISPENSER GUN  
CARTRIDGES

\*Source: Cluster Analysis - buying behavior and profiling, Keystone 2024  
\*\* Dynamic mix tips will be used with a bayonet ring



# A SOLUTION FOR EVERY SMILE.

Sultan Healthcare offers clinicians reliable solutions for the full spectrum of restoration needs—from impression to temporary crown and bridge materials. It's a synergistic system that delivers exceptional experiences and successful outcomes.



## Silgimix® alginate replacement impression material

This material combines the low cost of an alginate with the detail and versatility of a VPS. Silgimix® material works hand-in-hand with Versa-Temp® 2 material.



## Versa-Temp® temporary crown and bridge resin

From single-unit crowns to more demanding cases, Versa-Temp® and Versa-Temp® 2 enable beautiful temporaries that contribute to improved clinical results.



## SensiTemp® temporary cements

Available in an automix delivery system, SensiTemp® offers no mess, easy cleanup, minimal waste and maximum patient comfort.

## GENIE® PUTTY

Code	Setting Time	Pack
77655	Standard	300 ml Base + 300 ml Catalyst + 2 dosing spoons
77650	Fast	300 ml Base + 300 ml Catalyst + 2 dosing spoons

## GENIE® ORIGINAL 50 ml

Code	Viscosity	Setting	Pack
77635	Heavy Body	Standard	2 x 50 ml cartridge (Base + Catalyst) + 6 green mixing tips
77630	Heavy Body	Rapid	2 x 50 ml cartridge (Base + Catalyst) + 6 green mixing tips
77625	Regular Body	Standard	2 x 50 ml cartridge (Base + Catalyst) + 6 green mixing tips
77620	Regular Body	Rapid	2 x 50 ml cartridge (Base + Catalyst) + 6 green mixing tips
77615	Light Body	Standard	2 x 50 ml cartridge (Base + Catalyst) + 6 yellow mixing tips + 6 intraoral tips
77610	Light Body	Rapid	2 x 50 ml cartridge (Base + Catalyst) + 6 yellow mixing tips + 6 intraoral tips
77605	Extra Light Body	Standard	2 x 50 ml cartridge (Base + Catalyst) + 6 yellow mixing tips + 6 intraoral tips
77600	Extra Light Body	Rapid	2 x 50 ml cartridge (Base + Catalyst) + 6 yellow mixing tips + 6 intraoral tips
77715	Heavy Body	Standard	Bulk: 60 x 50 ml cartridge (Base + Catalyst)
77710	Heavy Body	Rapid	Bulk: 60 x 50 ml cartridge (Base + Catalyst)
77705	Regular Body	Standard	Bulk: 60 x 50 ml cartridge (Base + Catalyst)
77700	Regular Body	Rapid	Bulk: 60 x 50 ml cartridge (Base + Catalyst)
77725	Light Body	Standard	Bulk: 60 x 50 ml cartridge (Base + Catalyst)
77720	Light Body	Rapid	Bulk: 60 x 50 ml cartridge (Base + Catalyst)

## GENIE MAGIC MIX® 380 ml

Code	Viscosity	Setting	Pack
78835	Heavy Body	Standard	Bulk: 4 x 380 ml cartridge (Base + Catalyst) + 1 bayonet ring
78830	Heavy Body	Rapid	Bulk: 4 x 380 ml cartridge (Base + Catalyst) + 1 bayonet ring
78815	Heavy Body	Standard	1 x 380 ml cartridge (Base + Catalyst) + 10 dynamic mixing tips + 1 bayonet ring
78810	Heavy Body	Rapid	1 x 380 ml cartridge (Base + Catalyst) + 10 dynamic mixing tips + 1 bayonet ring
78805	Regular Body	Standard	1 x 380 ml cartridge (Base + Catalyst) + 10 dynamic mixing tips + 1 bayonet ring
78800	Regular Body	Rapid	1 x 380 ml cartridge (Base + Catalyst) + 10 dynamic mixing tips + 1 bayonet ring

## GENIE® BITE

Code	Setting Time	Pack
77640	Fast	2 x 50 ml cartridge (Base + Catalyst) + 6 green mixing tips
77645	Super Fast	2 x 50 ml cartridge (Base + Catalyst) + 6 green mixing tips

## SILGIMIX®

Code	Pack
79020	8 x 50 ml cartridge (Base + Catalyst)
79025	Bulk: 60 x 50 ml cartridge (Base + Catalyst)
79010	1 x 380 ml cartridge + 10 dynamic mixing tips + 1 bayonet ring
79015	Bulk: 4 x 380 ml cartridge + 1 bayonet ring

## IMPRESSION ACCESSORIES\*

Code	Pack
72001	48 x green mixing tips
72003	48 x yellow mixing tips
72007	100 x yellow intraoral tips
70051	Dispenser Gun (S-50) 50 ml
78840	50 x yellow Dynamic Mixing tips
77100	VPS Tray Adhesive

\* Refer to the website for the complete catalogue.

## VERSA-TEMP® 2

Code	Shade	Pack
72041	BL	1 x 50 ml cartridge + 15 light blue mixing tips
72051	A1	1 x 50 ml cartridge + 15 light blue mixing tips
72056	A2	1 x 50 ml cartridge + 15 light blue mixing tips
72061	A3	1 x 50 ml cartridge + 15 light blue mixing tips
72066	A1	Bulk: 3 x 50 ml cartridge
72071	A2	Bulk: 3 x 50 ml cartridge
72076	A3	Bulk: 3 x 50 ml cartridge

## VERSA-TEMP®

Code	Shade	Pack
72015	A1	1 x 25 ml cartridge + 10 light blue mixing tips
72020	A2	1 x 25 ml cartridge + 10 light blue mixing tips
72025	A 3.5	1 x 25 ml cartridge + 10 light blue mixing tips
72030	B1	1 x 25 ml cartridge + 10 light blue mixing tips
72050	A1	1 x 50 ml cartridge + 15 light blue mixing tips
72055	A2	1 x 50 ml cartridge + 15 light blue mixing tips
72060	A3	1 x 50 ml cartridge + 15 light blue mixing tips
72065	A1	Bulk: 3 x 50 ml cartridge
72070	A2	Bulk: 3 x 50 ml cartridge
72075	A3	Bulk: 3 x 50 ml cartridge

## SENSITEMP® RESIN

Code	Pack
70010	5 ml syringe + 10 brown mixing tips
70020	25 ml cartridge + 15 light blue mixing tips

## TEMP ACCESSORIES

Code	Pack
70025	Dispenser Gun (S-25) 25 ml
70051	Dispenser Gun (S-50) 50 ml
70035	50 x light blue mixing tips (SensiTemp® Resin 25 ml; Versa-Temp®; Versa-Temp® 2)
70030	20 x blunt brown mixing tips (SensiTemp® Resin 5 ml)



[www.sultanhealthcare.com](http://www.sultanhealthcare.com)